

INOGEN ONE G3 PRODUCT OVERVIEW

Provider/Distributor

oxygen. anytime. anywhere.

0

(0





- ✓ Setting 5 added = 1050 ml/min
- ✓ New noise reduction 39 dBA (on setting 2)
- Inlet vents are now above exhaust vents
- ✓ New, more rigid Gross Particle Filter
- New, higher capacity Zeolite, requiring differentiated columns: orange dust caps, blue sticker on end caps, new serial number format.





Inogen One G3 Concentrator

- One Single Battery (other battery configurations available)
- Nasal Cannula
- Carry Bag
- AC Power Supply
- DC Power Cable
- User Manual





Dimensions:	L/W/H: 8.75in./3.0in./7.25in	
With single battery:	L/W/H: 8.75in./3.0in./8.25in.	
Weight:	4.8 pounds (2.19Kg.) (includes single battery)	
Noise:	39 dBA (on setting 2)	
Warm-Up Time:	2 minutes	
Oxygen Concentration:	90% -3% / +6% at all settings	
Flow Control Settings:	5 settings: 1 to 5	
Power: AC Power Supply:	AC Input: 100 to 240VAC 50 to 60Hz	
DC Power Cable: Rechargeable Battery:	Auto-Sensing: 1.0A DC Input: 13.5-15.5VDC,10A Max. Voltage: 12.0 to 16.8VDC	
Battery Duration:	Up to 4.5 hours with single battery Up to 9 hours with double battery	
Battery Charging Time:	Up to 4 hours for a single battery Up to 8 hours for a double battery	
Environmental Ranges Intended for Use:	Temperature:41 to 104°F (5 to 40°C)Humidity:0% to 95%, non-condensingAltitude:0 to 10,000 ft (0 to 3048 meters)	



Battery Options





Model # BA-300 8 Cell Standard Battery ~ 4 Hours Run Time Included with System

Model # BA-316 16 Cell Optional Battery ~ 8 Hours Run Time Sold as Accessory





Battery Run Time

	Power Consumption	8 cell run time	16 cell run time
Setting 1.0	19 W	4.5 hours	9.5 hours
Setting 2.0	20 W	4.0 hours	8.5 hours
Setting 3.0	25 W	3.0 hours	6.0 hours
Setting 4.0	40 W	2.0 hours	4.0 hours
Setting 5.0	50 W	1.7 hours	3.3 hours

Optional Accessory External Battery Charger Model # BA-303



When the battery is in the correct position, a solid red light will indicate that the battery is charging.

When the green light illuminates, the battery is fully charged.



Battery Location







Latch will click into place when the battery is in position

Air Intake and Exhaust





Air Intake and Exhaust









Blue round stickers are applied to the end caps of the columns to differentiate. In addition, for a pair of replacement columns the dust caps will be orange color instead of black.

Access by removing battery and sliding latches to release columns.

CONTRACTOR OF BRIDE PROVIDE

Column serial numbers will include "HF" to differentiate as new columns.

Inogen One G3 Carry Bag Model # CA-300



G3 intake and exhaust vents should be visible through mesh













Inogen One G3 vs. SimplyGo Mini Inogen			
- Flow Settings	1 - 5	1 - 5	
Maximum Oxygen Output	1050 ml/min	1000 ml/min	
Sound Level @ Setting 2	39 dBA	43 dBA	
Size w/Single Battery	8.75 in. x 3.0 in. x 8.25 in.	9.4 in. x 3.6 in. x 8.3 in.	
Weight w/Single Battery	4.8 lbs.	5.0 lbs.	
Battery Charge Times	4 hours/8 hours	4 hours/8 hours	
Battery Duration	Up to 4.5 hours/9 hours	Up to 4.5 hours/9 hours	
FAA Approval	Yes	No	
oxygen. anytime. anywhere.			



Low Priority Alerts

The following low priority alert messages are accompanied by a **double beep** and a **solid red light**.

Current Message Display and Text	New G3 Message Display and Text	Condition/Action/Explanation
O2 Service Soon	Replace Columns	Column maintenance is required within 30 days. Contact your equipment provider to arrange for service.
Oxygen Low See Manual	Oxygen Low	Concentrator is producing oxygen at a slightly low level (<82%) for a period of 10 minutes. If condition persists, contact your equipment provider.
Service Needed	Service Soon	The concentrator requires servicing at the earliest convenience. The concentrator is operating to specification and may continue to be used. Contact your equipment provider to arrange for service.
O2 Sensor Fail See Manual	Sensor Fail	The concentrator's oxygen sensor has malfunctioned. You may continue to use the concentrator. If the condition persists, contact your equipment provider.

Inogen One G3 Error Code Changes Inogen

Medium Priority Alerts

The following low priority alert messages are accompanied by a **triple beep**, **repeated every 25 seconds**, and a **flashing red light**.

Current Message Display and Text	New G3 Message Display and Text	Condition/Action/Explanation
Oxygen Error Service Needed	Oxygen Error	Oxygen output concentration has been below 50% for 10 minutes. If condition persists, switch to your backup oxygen source and contact your equipment provider to arrange for service.
System Error See Manual	No message will display this code has been obsoleted as a medium priority and is replacing the High Priority text for Error (###)	Concentrator is producing oxygen at a slightly low level (<82%) for a period of 10 minutes. If condition persists, contact your equipment provider.

nogen One G3 Error Code Changes Inogen

High Priority Alerts

oxyger

The following high priority alert messages are accompanied by a five beep pattern, repeated every 10 seconds and a flashing red light.

	Current Message Display and Text	New G3 Message Display and Text	Condition/Action/Explanation
	Battery HOT Shut Down	Battery HOT	 Battery has exceeded temperature limit while concentrator is running on battery power. Concentrator has stopped producing oxygen. If possible, move concentrator to a cooler location, then turn power off and back on. Ensure air intake and outlet vents have clear access and particle filters are clean. If condition persists, switch to a backup source of oxygen and contact your equipment provider.
1111	System HOT Shut Down	System HOT	Concentrator temperature is too high and oxygen production is shutting down. Ensure air intake and outlet vents have clear access and particle filters are clean. If condition persists, switch to a backup source of oxygen and contact your equipment provider.
	System COLD Shut Down	System COLD	This may result from the concentrator being stored in a cold environment (below 0°C (32°F)). Move to a warmer environment to allow the unit to warm up before starting it. If condition persists, switch to a backup source of oxygen and contact your equipment provider.
	Error (###) Service Needed	System Error	Concentrator has stopped producing oxygen and is shutting down. You should: 1. Switch to backup oxygen source 2. Contact your equipment provider

Inogen One® G3 Software Changes



G3 Software upgrades were implemented in 2015 that were retroactively applied to all existing G3 units and will still apply to the new Inogen One® G3. The patient-facing changes that came with those updates are listed again here:

- The Audible Alarm Button will toggle the Inogen One® G3's breath detection audible alert on and off. If power is lost, the breath detection audible alert **remains set in the user preferred mode.**
- The Display Backlight Button, after pushed, automatically turns off after 15 seconds.
- The Replace Columns Low Priority Alert can no longer be cleared temporarily by the patient. The error will now alternately display
 "Replace Columns" and the default display every 30 seconds until the columns are replaced and reset.





INOGEN ONE G3 OXYGEN CONCENTRATOR

TECHNICAL MANUAL

96-03996-00-01

Revision F

Frequently Asked Questions



- Q: Can the New Sieve Beds be used in the 4 Settings G3 device?
- A: No
- Q: Will the New Sieve Beds last longer with the Higher Capacity Zeolite?
- A: No
- Q: Is the Carrying Bag the same or is there a new carrying bag for the G3 HF?
- A: The bag has not changed.
- Q: Is the weight the same or slightly lower due to the internal changes and New Zeolite?
- A: The weight is slightly less averaging at 4.76 with single battery.
- Q: Are the Air Inlets and Outlets bigger?
- A: Yes they are slightly larger to allow more air in.



THANK YOU!