

Neb-u-Tyke Train Nebulizer Compressor

- Reliable, Piston-driven Compressor
- Colorful, Child-Friendly design
- Interactive stickers for patient to customize their train
- Comes complete with Nebulizer Accessory Kit
- 5 year limited warranty
- Designed to meet the requirements of HCPCS Code: E0570

JB0112-164 Neb-u-Tyke Train
JB77300 Neb-u-Tyke Train & Depot Case
JB76230 Depot Case
TRAIN-DISPLAY Display



The John Bunn Neb-u-Tyke Train Nebulizer comes complete with:

(single patient use only)

- 1 adult mask
- 1 pediatric mask
- 1 mouthpiece
- 1 PVC tubing
- 1 "T" adapter
- 3 spare filters
- 1 reservoir tube



Lights up during use!



Train Display

Technical Specifications

Power	115VAC, 60Hz
Power Consumption	≤ 90 watts
Operating Temperature Range	41°F TO 104°F
Operating Humidity Range	30 to 95% RH
Storage Temperature Range	14°F TO 158°F
Storage Humidity Range	10 to 95% RH
Dimensions	9.65" (L) × 5.5" (H) × 4.5" (W)
Weight	3.5 lbs
Max Pressure (Compressor)	32.5 psi
Max Flow (Compressor)	Approx 10 LPM
Particle Size (Nebulizer)	≤5 microns
Capacity (Nebulizer)	5 ml
Method Of Nebulization	Piston Driven
Sound Level	under 55 dBA

Specifications are subject to change without notice.



GRAHAM-FIELD